


## Product Datasheet

### TXN Human Over-expression Lysate (orb1365796)

|                        |   |
|------------------------|---|
| <b>Description</b>     | TXN HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | TXN   |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>H9ZYJ2, P10599</b>   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |



Western blot validation of Left: Cell ly...