

## Product Datasheet

# TWEAK (TNFSF12) Human Over-expression Lysate (orb1366056)

|                        |   |
|------------------------|---|
| <b>Description</b>     | TNFSF12 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | TNFSF12   |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <a href="#">Q4ACW9</a> , <a href="#">O43508</a>                               |
| <b>Expiration Date</b> | 6 months from date of receipt.  |

