


## Product Datasheet

### FIBP Human Over-expression Lysate (orb1366259)

|                        |  |
|------------------------|--|
| <b>Description</b>     | FIBP HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | FIBP   |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">O43427</a>   |
| <b>Expiration Date</b> | 6 months from date of receipt.   |



Western blot validation of Left: Cell ly...