


## Product Datasheet

# DGCR2 Human Over-expression Lysate (orb1366705)

|                        |   |
|------------------------|---|
| <b>Description</b>     | DGCR2 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | DGCR2   |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>Q8IWC8, P98153</b>   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |



Western blot validation of Left: Cell ly...