



Product Datasheet

ABCC4 antibody (orb136696)


Description	Rabbit polyclonal antibody to ABCC4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human MRP4/ABCC4(NP_005836.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140-200KDa
Uniprot ID	O15439
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of extracts of var...



IHC-P staining of human esophageal cancer...