


Product Datasheet

Tubulin alpha antibody (orb137121)

Description	Mouse monoclonal to Tubulin alpha. Tubulin is one of several members of a small family of globular p...
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB:1:3000:1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of skeletal muscle...