



Product Datasheet

NPTX1 antibody (orb137172)

Description	Goat polyclonal to NPTX1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	neuronal pentraxin 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human Brain (Amygdala and Frontal Cortex) lysates (calculated MW of 47.1kDa according to NP_002513.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour. Additional bands of unknown identity was also consist. IHC: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 5µg/ml.
Protein Sequence	ERVKIETALTS LHQR
Clonality	Polyclonal
MW	47.1
Entrez	4884
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Frontal C...



Immunohistochemical staining of Human Ce...