



## Product Datasheet

### HOXA5 antibody (orb137175)

<b>Description</b>	Goat polyclonal to HOXA5
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Bovine, Human, Mouse, Porcine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	HOXA5 (aa83-92)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: Approx 28kDa band observed in Human, Mouse, Rat and Pig Kidney lysates (calculated MW of 29.2kDa according to Mouse NP_034583.1). Recommended concentration: 1-3µg/ml.
<b>Protein Sequence</b>	EPRYSQPATS
<b>Clonality</b>	Polyclonal
<b>MW</b>	29.3
<b>Entrez</b>	<a href="#">3202</a>
<b>Dilution Range</b>	ELISA: 1:8000, WB: 1-3 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of staining of Hum...



Western blot analysis of staining of Pig...