


## Product Datasheet

# Elongation of very long chain fatty acids protein 3 (ELOVL3) Human Over-expression Lysate (orb1373634)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | ELOVL3 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | ELOVL3                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>Q9HB03</b>                                                                |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western blot validation of Left: Cell ly...