

## Product Datasheet

### ICAM-3 antibody (orb137661)

**Description**

Mouse monoclonal to ICAM-3. Intercellular

CD50 (Cluster of...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FC, ICC, IHC-Fr, IP, WB

**Preservatives**

Each vial contains 100 µl 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

**Form/Appearance**

Each vial contains 100 µl 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

**Storage**

Storage: The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product

**Note**

For research use only

**Isotype**

IgG1/k

**Clonality**

Monoclonal

**Source**

BH10 is a mouse monoclonal IgG1/k antibody derived by fusion of mouse Ag8.653 cells with spleen cells from a BALB/c mouse immunized with an ICAM-3/HEK transfectant.

**Uniprot ID**

[P32942](#)

**Expiration Date**

12 months from date of receipt.