


Product Datasheet

Human CA9/Carbonic anhydrase 9 ELISA Kit (orb138063)

Description	Human CA9/Carbonic anhydrase 9 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: P59-D414
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-259 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	49698 MW
Uniprot ID	Q16790
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CA9
PicoKine
ELISA Kit
standard
cu...