


Product Datasheet

Human ANGPTL4 / Angiopoietin Like 4 / ARP4 ELISA Kit (orb138079)

Description	Human ANGPTL4 / Angiopoietin Like 4 / ARP4 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L165-S406
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 1.2 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	45214 MW
Uniprot ID	Q9BY76
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human ANGPTL4 PicoKine ELISA Kit standar...