


Product Datasheet

Human MICB ELISA Kit (orb138082)

Description	Human MICB ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-G298
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	42646 MW
Uniprot ID	Q29980
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human MICB
PicoKine
ELISA Kit
standard c...