


Product Datasheet

Human IL-17RC ELISA Kit (orb138088)

Description	Human IL-17RC ELISA Kit
Reactivity	Human
Range	250 pg/ml - 16000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L21-A454
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	86240 MW
Uniprot ID	Q8NAC3
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Human IL-17RC PicoKine ELISA Kit standard...