


Product Datasheet

Mouse Complement C5a ELISA Kit (orb138090)

Description	Mouse Complement C5a ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N679-R755
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is around 10 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	188878 MW
Uniprot ID	P06684
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Complement C5a PicoKine ELISA Kit ...