


Product Datasheet

Human IGFBP7/Igfbp Rp1 ELISA Kit (orb138093)

Description	Human IGFBP7/Igfbp Rp1 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D30-L282
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 50 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	29130 MW
Uniprot ID	Q16270
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
IGFBP7
PicoKine
ELISA Kit
standard...