


Product Datasheet

Human Galectin-9 / Gal 9 ELISA Kit (orb138106)

Description	Human Galectin-9 / Gal 9 ELISA Kit
Reactivity	Human
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-T355
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:5, concentration in serum and plasma is 2.48-12.48 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39518 MW
Uniprot ID	000182
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Galectin-9 PicoKine ELISA Kit stan...