

Product Datasheet

MIP1 gamma antibody (FITC) (orb13834)

Description	MIP1 gamma antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from mouse CCL9
Target	Macrophage Inflammatory Protein 1 gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5FVN3
Expiration Date	12 months from date of receipt.