

Product Datasheet

Clenbuterol Hydrochloride antibody (orb1394)

Description	Clenbuterol Hydrochloride antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Clenbuterol
Target	Clenbuterol
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	0.31365 kDa
CAS Number	21898-19-1
Expiration Date	12 months from date of receipt.