



## Product Datasheet

### THAP11 antibody (orb140008)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse Monoclonal to THAP11.   |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | Synthetic Peptide   |
| <b>Target</b>              | THAP11  |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                            |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Gene Symbol: THAP11. Gene-IDs: 57215. Accessions: Q96EK4, NP_065190.2.  |
| <b>Isotype</b>             | IgM   |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 34455   |
| <b>Uniprot ID</b>          | <a href="#">Q96EK4</a>  |
| <b>NCBI</b>                | <a href="#">NP_065190.2</a>   |
| <b>Dilution Range</b>      | WB: 1:100-500   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot analysis of THAP11 antibody...