



## Product Datasheet

### PLS1 antibody (orb140065)

<b>Description</b>	Rabbit Polyclonal to PLS1.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	PLS1
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Gene Symbol: PLS1. Gene-IDs: 5357. Accessions: Q14651, NP_001165783.1, NP_001138791.1.
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	70253
<b>Uniprot ID</b>	<a href="#">Q14651</a>
<b>NCBI</b>	<a href="#">NP_002661.2</a> , <a href="#">NP_001165783.1</a> , <a href="#">NP_001138791.1</a>
<b>Dilution Range</b>	WB: 1:1000, IHC-P: 1:10-50
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical analysis of formalin...