



Product Datasheet

DYSFIP1 antibody (orb140075)

Description	Rabbit Polyclonal to DYSFIP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PPP1R27
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: DYSFIP1. Gene-IDs: 116729. Accessions: NP_001007534.1, Q86WC6.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17438
Uniprot ID	Q86WC6
NCBI	NP_001007534.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...