



Product Datasheet

MRAP2 antibody (orb140084)

Description	Rabbit Polyclonal to MRAP2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	MRAP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: MRAP2. Gene-IDs: 112609. Accessions: NP_612418.2, Q96G30.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23548
Uniprot ID	Q96G30
NCBI	NP_612418.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...