



Product Datasheet

PGAP3 antibody (orb140103)

| | |
|----------------------------|--|
| Description | Rabbit Polyclonal to PGAP3. |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Immunogen | Synthetic Peptide |
| Target | PGAP3 |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Application notes | Gene Symbol: PGAP3. Gene-IDs: 93210. Accessions: Q96FM1, NP_219487.3. |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 36475 |
| Uniprot ID | Q96FM1 |
| NCBI | NP_219487.3 |
| Dilution Range | WB: 1:1000, IHC-P: 1:10-50 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of formalin...