



Product Datasheet

ASB4 antibody (orb140106)

Description	Rabbit Polyclonal to ASB4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ASB4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: ASB4. Gene-IDs: 51666. Accessions: Q9Y574, NP_057200.1, NP_665879.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48195
Uniprot ID	Q9Y574
NCBI	NP_057200.1 , NP_665879.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...