



Product Datasheet

CHMP2B antibody (orb140120)

Description	Rabbit Polyclonal to CHMP2B.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CHMP2B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: CHMP2B. Gene-IDs: 25978. Accessions: Q9UQN3, NP_054762.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23907
Uniprot ID	Q9UQN3
NCBI	NP_054762.2 , NP_001231573.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...