



Product Datasheet

MOSPD2 antibody (orb140127)

Description	Rabbit Polyclonal to MOSPD2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	MOSPD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: MOSPD2. Gene-IDs: 158747. Accessions: NP_001170946.1, NP_689794.1, Q8NHP6.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59746
Uniprot ID	Q8NHP6
NCBI	NP_689794.1 , NP_001170946.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...