



Product Datasheet

Mouse TLR6 antibody (orb140161)

Description	Rabbit Polyclonal to Mouse TLR6.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	Tlr6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: Tlr6. Gene-IDs: 21899. Accessions: NP_035734.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91116
Uniprot ID	Q9EPW9
NCBI	NP_035734.3
Dilution Range	IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse TLR6 anti...