

Product Datasheet

Anti-ChAT (Choline Acetyltransferase) Antibody (orb1401953)

Description	Our Anti-ChAT (Choline Acetyltransferase) primary antibody from Biorbyt is goat polyclonal. It detects chicken, guinea pig, human, mouse, and non-human primate, and rat ChAT (Choline Acetyltransferase) and is affinity purified.
Species/Host	Goat
Reactivity	Gallus, Guinea pig, Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant human choline acetyltransferase
Target	ChAT (Choline Acetyltransferase)
Preservatives	PBS + 5 mg/ml rAlbumin + 0.2% NaN ₃
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
MW	70
Uniprot ID	P28329

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Dilution Range	WB: 1:100, IHC: 1:100-1:1000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)