

Product Datasheet

NR5A2 antibody (orb140290)

Description	Rabbit Polyclonal to NR5A2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NR5A2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: NR5A2. Gene-IDs: 2494. Accessions: NP_003813.1, O00482, NP_995582.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61331
Uniprot ID	O00482
NCBI	NP_001263393.1 , NP_995582.1 , NP_003813.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

