

Product Datasheet

SAA4 antibody (orb140347)

Description	Rabbit Polyclonal to SAA4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SAA4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: SAA4. Gene-IDs: 6291. Accessions: NP_006503.1, P35542.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14747
Uniprot ID	P35542
NCBI	NP_006503.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

