

Product Datasheet

HNRNPH3 antibody (orb140371)

Description	Rabbit Polyclonal to HNRNPH3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HNRNPH3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: HNRNPH3. Gene-IDs: 3189. Accessions: NP_036339.1, P31942, NP_067676.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36926
Uniprot ID	P31942
NCBI	NP_067676.2 , NP_036339.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

