

Product Datasheet

NDRG4 antibody (orb140379)

Description	Rabbit Polyclonal to NDRG4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NDRG4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: NDRG4. Gene-IDs: 65009. Accessions: Q9ULP0, NP_001123959.1, NP_065198.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38459
Uniprot ID	Q9ULP0
NCBI	NP_001229762.1 , NP_001229765.1 , NP_065198.1 , NP_001123959.1 , NP_001229764.1 , NP_001229763.1 , NP_075061.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

