

Product Datasheet

RPL10A antibody (orb140385)

Description	Rabbit Polyclonal to RPL10A.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RPL10A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: RPL10A. Gene-IDs: 4736. Accessions: NP_009035.3, P62906.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24831
Uniprot ID	P62906
NCBI	NP_009035.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

