

Product Datasheet

YAP (Phospho-Ser127) antibody (orb140432)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal to YAP (Phospho-Ser127) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | DOT |
| Immunogen | Synthetic Peptide |
| Target | YAP1 |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Application notes | Gene Symbol: YAP1. Gene-IDs: 10413. Accessions: P46937, NP_001123617.1. |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 54462 |
| Uniprot ID | P46937 |
| NCBI | NP_001269029.1 , NP_001269030.1 , NP_001269027.1 , NP_001181973.1 , NP_001181974.1 , NP_001269028.1 , NP_001123617.1 |
| Dilution Range | DOT: 1:500 |
| Expiration Date | 12 months from date of receipt. |

