

## Product Datasheet

### PAFAH1B2 (Phospho-Ser64) antibody (orb140435)

<b>Description</b>	Rabbit polyclonal to PAFAH1B2 (Phospho-Ser64)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	DOT
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	PAFAH1B2
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Gene Symbol: PAFAH1B2. Gene-IDs: 5049. Accessions: NP_001171677.1, NP_002563.1, P68402.
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	25569
<b>Uniprot ID</b>	<b>P68402</b>
<b>NCBI</b>	<b>NP_001171675.1, NP_002563.1, NP_001171677.1, NP_001171676.1</b>
<b>Dilution Range</b>	DOT: 1:500
<b>Expiration Date</b>	12 months from date of receipt.

