

Product Datasheet

MYH9 (Phospho-Tyr158) antibody (orb140437)

Description	Rabbit polyclonal to MYH9 (Phospho-Tyr158)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MYH9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	Gene Symbol: MYH9. Gene-IDs: 4627. Accessions: P35579, NP_002464.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	226532
Uniprot ID	P35579
NCBI	NP_002464.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

