
Product Datasheet

**Penicillin (Group of beta-Lactam Antibiotics)
antibody (orb1410468)**

Description

Mouse monoclonal antibody to

Antibiotics)

Species/Host

Mouse

Conjugation

Unconjugated

Tested

ELISA

Applications**Immunogen**

Penicilloyl-transferrin conjugate

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide store at 2 to 8 °C. Antibody without azide store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

356 Daltons

Uniprot ID

Not Applicable

Entrez

Not Applicable

Dilution Range

ELISA (For coating, use Ab at 1-5ug/ml; order Ab without BSA)Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.