

Product Datasheet

BTN3A1/2/3 antibody (orb1410829)

Description	Rabbit polyclonal antibody to BTN3A1/2/3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O00481 , P78410 , O00478
Dilution Range	IHC: IHC-p 1:50-200, ELISA 1:10000-20000
Expiration Date	12 months from date of receipt.