

## Product Datasheet

### TLE1/2/3/4 antibody (orb1411027)

<b>Description</b>	Rabbit polyclonal antibody to TLE1/2/3/4.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q04727</a> , <a href="#">Q04726</a> , <a href="#">Q04724</a> , <a href="#">Q04725</a>
<b>Dilution Range</b>	WB: WB 1:500-2000, ELISA 1:10000-20000
<b>Expiration Date</b>	12 months from date of receipt.