

Product Datasheet

Sox-8/9/17/18 antibody (orb1412129)

Description	Rabbit polyclonal antibody to Sox-8/9/17/18.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P48436, P57073, Q9H6I2, P35713
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.