

## Product Datasheet

### Olfactory receptor 10X1 antibody (orb1412782)

<b>Description</b>	Rabbit polyclonal antibody to Olfactory receptor 10X1.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q8NGY0</a> , <a href="#">Q8NGX7</a>
<b>Dilution Range</b>	WB: Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.