

## Product Datasheet

### Galectin-9 antibody (orb1413657)

|                            |                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Galectin-9.                                                                                                                              |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                           |
| <b>Tested Applications</b> | IF, IHC-P, WB                                                                                                                                                          |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                 |
| <b>Note</b>                | For research use only                                                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                             |
| <b>Uniprot ID</b>          | <b>O00182</b>                                                                                                                                                          |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>Immunohistochemistry: 1/100 - 1/300.<br>Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                        |