

Product Datasheet

Na⁺/K⁺-ATPase α 1 (Phospho-Tyr260) antibody (orb1414835)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase α 1 (Phospho-Tyr260).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05023
Dilution Range	WB: WB 1:500-2000, ELISA 1:10000-20000
Expiration Date	12 months from date of receipt.