

Product Datasheet

Histone H2A (phospho Thr121) antibody (orb1414889)

Description	Rabbit polyclonal antibody to Histone H2A (phospho Thr121).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96KK5 , Q7L7L0 , P20671 , P04908 , Q6FI13 , P0C0S8
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.