

Product Datasheet

p70 S6 kinase α (phospho Ser427) antibody (orb1414894)

Description	Rabbit polyclonal antibody to p70 S6 kinase α (phospho Ser427).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P23443
Dilution Range	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.