

## Product Datasheet

### Trk A (phospho Tyr791) antibody (orb1414918)

<b>Description</b>	Rabbit polyclonal antibody to Trk A (phospho Tyr791).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P04629</b>
<b>Dilution Range</b>	IHC: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.