

Product Datasheet

CD71 (phospho Ser24) antibody (orb1415029)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to CD71 (phospho Ser24). |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P02786 |
| Dilution Range | WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. |
| Expiration Date | 12 months from date of receipt. |