

Product Datasheet

Arrestin- β -1 (phospho Ser412) antibody (orb1415073)

Description	Rabbit polyclonal antibody to Arrestin- β -1 (phospho Ser412).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P49407
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.