

## Product Datasheet

### EphA2/5 (phospho Tyr594) antibody (orb1415173)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to EphA2/5 (phospho Tyr594).  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">P54756</a> , <a href="#">P29317</a>  |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>ELISA: 1/10000. Not yet tested in other applications.   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |